
Meisinger

since
1888

PODIATRY & FOOT CARE CATALOGUE

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The world is evolving faster and faster all the time. Do you get that feeling, too? What was new yesterday is old hat today. Innovation is a central task for us also as a manufacturer of precision instruments “Made in Germany“. Always one step ahead – that’s our distinguishing feature since 1888.

Hager & Meisinger GmbH is one of the worldwide leading developers and manufacturers of rotating high-tech instruments in the field of foot care, podiatry and dental medicine. The overall assortment includes over 12,000 products. At our headquarters located in Neuss, around 45 million individual instruments are produced annually, which are sold in more than 100 countries under the brand MEISINGER.

Our portfolio in the field of rotating instruments is consistently aligned with the requirements of daily practice in the field of foot care and podiatry. Effective and gentle removal with a high level of instrument durability – that is your requirement, and constitutes our motivating force.

The new edition of our podology catalogue fills us with pride – instruments that have not only proven their worth many times over but also the newest “state of the art” instruments that Meisinger strives to provide worldwide.

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DIAMOND INSTRUMENTS (SHANK 104)

MEISINGER diamond instruments basically consist of high-strength, stainless-steel basic profiles and shanks. By means of the most modern coating process only natural and selected diamond grains are used. The result is a uniform, safe and effective diamond coating in five different grain sizes.

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
Grits of the Diamond Instruments and relevant area of application

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

DIAMOND INSTRUMENTS (SHANK 314)

MEISINGER diamond instruments with FG shanks (Ø 1,6 mm) are designed for the usage in turbines. The water cooling in combination with high speeds enables efficient work and long service life of the instruments. Please note that instruments with FG shanks can only be used in appropriate turbines. The use in conventional handpieces is not possible.

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

TUNGSTEN CARBIDE CUTTERS

MEISINGER tungsten carbide cutters are produced from highly selected compressed fine grain tungsten carbide. Its firm construction and highly precise toothing ensures optimal accuracy to size and concentricity.

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Toothing carbide metal cutters and corresponding fields of application

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Double-gearing, medium (DG) and fine (DF)

Shank HP		
Ref.-No.	HM251DG	HM251DF
Ø mm	060	060
	2	2

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Thanks to an innovative and completely novel cutting system, the time of making compromises is over: The new cutter with a twin blade combines the durability typical of carbide metal with efficient abrasiveness. The secret is a twin blade with blades moving at various depths. A system that many know from beard shaving and that has been tried and tested within the scope of such applications.

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The Callus Cleaner for all cases: The new Carbide Cutter with file-toothing ensure efficient removal of all types of callus skin, leaving absolutely smooth skin. It is no longer necessary to use different instruments depending on the condition and thickness of the cornea. This saves time and money!

The new Callus Cleaner are effective gadgets for everyday use and convince with high-quality working results. A must have for the modern foot care!

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

STEEL INSTRUMENTS (SHANK 104)

MEISINGER steel instruments are selected from stainless-steel and produced in one piece. Its optimal firm construction and its expedient, above average concentricity, precision tothing and dimensions make it a convincing instrument. It offers best possible material hardness combined with elasticity and is therefore ideal for use in podiatry.

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Type of Cuts and relevant area of application

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

82RF & 92 RF

- Specially for Mycose-nails
- Fast and gentle operation
- For thinning the nail-plate

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TREPHINES


Meisinger Trephines are manufactured from selected stainless steel. Its above average concentricity, precise, sharp cutting edges and its appropriate and stable construction make it a convincing instrument. Optimal hard material allows for careful removal of painful deep corns as well as removing and cutting out hard callus skin.

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
CERAMIC ABRASIVES

MEISINGER abrasives are particularly characterized by uniform toughness and grit, as well as by high strength, stainless steel and an exact measured shank. These properties guarantee a precise concentricity.

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Abrasive, pink (for fine work on corns and nails as for work on orthosis)

Shank HP														
Ref.-No.	601	602	603	613	618	620	623	624	626	627	637	638	639	640
Ø mm	030	040	050	070	050	050	060	060	050	060	020	025	035	050
L mm	3,0	4,0	5,0	3,0	1,0	3,0	3,0	5,0	1,5	1,5	6,0	6,0	6,0	6,0
	5	5	5	5	5	5	5	5	5	5	5	5	5	5

Shank HP														
Ref.-No.	648	649	650	652R	661	662	663	665	666	667	671	731	732	733
Ø mm	020	025	028	035	025	035	060	060	025	035	060	065	050	035
L mm	6,0	6,0	6,0	10,0	7,0	7,5	10,5	12,0	6,5	7,0	10,0	12,5	12,5	10,5
	5	5	5	5	5	5	5	5	5	5	5	5	5	5

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Abrasive, pink (for wide area grinding of callus skin)

Shank HP			
Ref.-No.	743S	744	751
Ø mm	130	110	100
L mm	22,0	17,0	20,0
	5	5	5

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Abrasive, white (for removal and smoothing of callus skin)

Shank HP	
Ref.-No.	453
Ø mm	100
L mm	10,0
	5

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POLISHERS

MEISINGER polishers consist of high quality polishing agents in a resistant, flexible connection of polysiloxan and indian rubber which are mounted on stainless steel shanks. These qualities guarantee a clean, orderless and shock-free polish by low heat development.

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



1435 Kit of Carbide Chip Coated for callus skin

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Observe the regulations precisely of the disinfectant manufacturer for aluminum products: some are not suitable!

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Stainless steel Bur Blocks, non-corrosive

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Designed for fast and easy cleaning of diamond instruments, for improved performance (Cleaning stone to be moistened thoroughly before use).

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Rust preventing spray, sterilizable, for steel instruments. The RS-Spray prevents rusting during and after sterilisation. Applicable for all methods of sterilisation.

Autoclave: max. 134°C / 273.2 °F; Hot air sterilisation: max. 180 °C / 356 °F; Contents: 200 ml (219 g)

	grob / coarse
Ref.-No.	529 MG
Dimensions	75 x 25 x 15 mm

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Rotating instruments by MEISINGER cover the entire spectrum of applications in podiatry. The respective maximum speed is a theoretical value given on the basis of the geometry of the instruments and the principal application. You can find the specific maximum and recommended speeds in the product-related information.

To over exceed the given rotation speed can lead to impairment of the efficiency of the products and/or be of danger for the patient as well as the user. The choice of the working speed for a specific application case is determined by the material subsequently to be used, the instrument power, and the maximal down force pressure ultimately decided upon by the user.

Podologie / Fußpflege Podiatry / Foot care		Podologie / Fußpflege Podiatry / Foot care		Podologie / Fußpflege Podiatry / Foot care	
Diamant diamond		Hartmetall / Diamant Tungsten carbide / diamond		Keramische Scheifer Abrasives	
Größe Size	FG FG	Größe Size	Handstück HP	Größe Size	Handstück HP
min ⁻¹ (M/S)		min ⁻¹ (M/S)		min ⁻¹ (M/S)	
005	300.000	010	50.000	025	50.000
006	300.000	005	50.000	027	50.000
007	300.000	006	50.000	029	50.000
008	300.000	007	50.000	031	50.000
009	300.000	008	50.000	033	50.000
010	300.000	009	50.000	035	50.000
012	300.000	010	50.000	037	50.000
014	300.000	012	50.000	040	50.000
016	280.000	014	50.000	042	50.000
018	250.000	016	50.000	045	50.000
021	210.000	018	50.000	047	50.000
023	190.000	021	50.000	050	50.000
025	180.000	023	50.000	055	50.000
027	160.000	025	50.000	060	50.000
029	150.000	027	50.000	065	50.000
031	150.000	029	50.000	070	50.000
033	120.000	031	50.000	075	50.000
035	120.000	033	50.000	080	50.000
037	120.000	035	50.000	085	50.000
040	100.000	037	50.000	090	45.000
042	100.000	040	50.000	095	45.000
		042	50.000	100	40.000
		045	50.000	110	35.000
		047	45.000	120	35.000
		050	40.000	130	30.000
		055	35.000		
		060	35.000		
		065	30.000		
		070	30.000		
		075	25.000		
		080	25.000		
		085	25.000		
		090	25.000		
		095	20.000		
		100	20.000		
		110	20.000		
		120	18.000		
		130	18.000		
		140	16.000		
		150	16.000		
		160	14.000		
		170	14.000		
		180	12.000		
		190	12.000		
		200	10.000		
		220	10.000		

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Please pay attention to the reference speed indicated on the packaging!

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Instructions for the reprocessing (cleaning, disinfecting and sterilizing) of medical devices produced by Hager & Meisinger GmbH Issued: August 2013

The medical devices produced and sold by Hager & Meisinger GmbH are reusable unless their label contains explicit information to the contrary. However, as a rule, it is the sole responsibility of the doctor/expert using the devices to decide whether, depending on the respective case and the potential wear and tear of the products, he can re-use the products and how frequently he uses them. In case of doubt, it is always advisable to discard the products early and to replace them. The manufacturer Hager & Meisinger GmbH cannot guarantee the faultless function and performance of the products combined with a maximum degree of safety if the products are overused.

These reprocessing instructions apply in principle to all medical devices making up the product range supplied by Hager & Meisinger GmbH. Any particular features and/or exclusions that only concern individual items or groups of items are referred to separately. As to the general application and safety instructions concerning the use of the products, we would advise you to consult the application and safety instructions for the medical products of Meisinger, which are available separately (see also at www.meisinger.de).

General rules

All products need to be cleaned, disinfected and sterilised prior to use; this applies in particular to the first-time use of products i.e. after their delivery because all products are supplied non-sterile (cleaning and disinfecting after the removal of transport packaging; sterilisation after removing wrapping). The thorough cleaning and disinfecting is indispensable in order to achieve effective sterilisation.

As you are responsible for the sterility of the products during use, please ensure

- that, as a rule, only sufficiently validated device and product-specific procedures are used for the cleaning/disinfecting and sterilisation process,
- that the used devices (disinfector, sterilizer) are maintained and checked at regular intervals and
- that the validated parameters are adhered to during each cycle.

Please ensure that, during use, you separate the soiled products and do not return them to the bur block/instrument tray to avoid further contamination of the loaded bur block/instrument tray. After cleaning/disinfecting the soiled products, arrange them in the bur block/instrument tray and sterilise the fully loaded bur block/instrument tray.

Please also take note of the statutory regulations valid in your country as well as the hygiene rules followed by the doctor's practice or the hospital. This applies in particular to the different rules regarding effective prion deactivation. You will find further instructions for the reprocessing of medical products and information on suitable products to carry out the reprocessing on the internet at www.dghm.de or at www.rki.de.

As the products are destined to be used for surgical, parodontological or endodontic procedures such as root canal debridement, they may penetrate the skin or the mucosa and come into contact with blood, internal tissues or organs (including wounds). Therefore, we recommend that they be assigned to risk group Critical B if used for their intended purpose.

Attention: For some products, additional aspects need to be considered (see chapter on 'Special instructions').

Cleaning and disinfecting

Basic rules:

Ideally, a mechanical procedure (disinfector) should be used for the cleaning and disinfecting process. Due to their clearly reduced effectiveness and reproducibility, manual procedures – even if they include the use of an ultrasound bath – should only be used if an automated procedure is unavailable. Pre-treatment needs to be carried out in both cases.

Apply hygiene-effective measures in line with the RKI guidelines at your workplace and carry out all work using powder and latex-free gloves.

Pre-treatment:

Abrasive impurities need to be removed from the products directly after use (within two hours maximum). To do so, use running water or a disinfectant solution; the disinfectant should be non-protein-fixing, its effectiveness should be established (e.g. it should be certified by the VAH (German Association for Applied Hygiene), the DGHM (German Society for Hygiene and Microbiology) or the FDA or have CE marking), it should be suitable for the disinfection of the products and compatible with the products (see chapter on 'Material resistance'). For the manual removal of impurities, only use nylon brushes intended for the purpose (also see the chapter on 'Special Instructions') or clean soft and lint-free cloths that you only use for this purpose. Do not use metal brushes or steel wool.

As to products with lumens (cavities): rinse all cavities three times by using a disposable syringe (minimum volume 5-10 ml) and a cannula.

Please note that the disinfectants used during pre-treatment only ensure personal protection and can be no substitute for the disinfection procedure to be used later - after completion of the cleaning process.

Mechanical cleaning/disinfection (disinfector/RDG (cleaning and disinfection device):

When choosing a disinfector you will have to ensure

- that the effectiveness of the disinfector has been certified (e.g. it has been licensed by the DGHM or the FDA or has CE marking according to DIN EN ISO 15883),
- that, if at all possible, a programme is used that has been certified for thermic disinfection (A0-value > 3000 or – with regard to older devices – at least 5 minutes at 90 °C (194 °F)) (chemical disinfection runs the risk of disinfectant residues remaining on the instrument),
- that the used programme is suitable for the products and has a sufficient number of rinsing cycles,
- that only sterile water or water with low levels of germs (max. 10 germs/ml) and endotoxins (max. 0.25 endotoxin units/ml) is used for the post-purge cycle (e.g. purified water/highly purified water),
- that the air used for drying is filtered and that the disinfector is regularly maintained and checked.

When choosing an appropriate cleaning system you need to ensure:

- that it is generally suitable for the cleaning of products made of metal and plastic,
- that, in addition, – if no thermic disinfection is used - a disinfectant is used whose efficacy has been certified (e.g. VAH/DGHM FDA licence or CE marking) and
- that it is compatible with the used cleaning agents and that the used chemicals are compatible with the products (see chapter on 'material resistance').

It is essential that the concentrations recommended by the manufacturer of the

cleaning and the disinfectant agent (if required) are adhered to at all times.

Procedure:

1. Disassemble the products as far as possible.
2. When using products with lumen: ensure that all lumens are rinsed effectively as part of the pre-cleaning process.
3. Place the disassembled products in the disinfector. Take care that the products do not touch. If possible: connect all instrument lumens to the rinsing attachment of the disinfector by using an appropriate rinsing adapter.
4. Start the programme.
5. Remove the products from the disinfector after the programme has finished.
6. Check and wrap the products straight after removal if possible (see chapters on 'Checking', 'Maintenance' and 'Packaging') if necessary, after they have been dried off completely in a clean place.

Proof of the general suitability for effective mechanical cleaning and disinfecting has been provided by an independent certified test laboratory using the disinfector G 7836 GD (thermic disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the cleaning agent Neodisher medizym (Dr. Weigert GmbH & Co. KG, Hamburg). As part of the testing the laboratory used the above described procedure.

Manual cleaning and disinfection:

When choosing an appropriate cleaning and disinfecting agent you need to ensure

- that they are generally suitable for the cleaning and/or disinfection of products made of metal and plastic,
- that the cleaning agent, where used, is suitable for ultrasound cleaning (no production of foam)
- that an appropriate disinfectant is used whose efficacy has been certified (e.g. VAH/DGHM FDA licence or CE marking) and that it is compatible with the used cleaning agent and
- that the used chemicals are compatible with the products (see chapter on 'Material resistance').

Ideally, combined cleaning/disinfecting agents should not be used. Combined cleaning/disinfecting agents can only be used in cases where there is a very low degree of contamination (no visible soiling).

It is essential that the concentrations and contact times recommended by the manufacturer of the cleaning and the disinfectant agents are adhered to at all times.

Only use freshly made solutions and water that is either sterile or low in germs (max. 10 germs/ml) and endotoxins (max. 0.25 endotoxin units/ml) (e.g. purified water/highly purified water) and only use filtered air for drying.

Procedure: Cleaning

1. Disassemble the products as far as possible.
2. Place the disassembled products in the cleaning bath for the prescribed contact time and ensure that products are sufficiently submerged; where required, use ultrasound support but at least careful brushing using the nylon brushes intended for this purpose (see also chapter on 'Special Instructions'). Take care that the products do not touch. With regard to products with lumens: rinse all instrument lumens five times at the beginning and or at the end of the contact time using a disposal syringe (minimum volume 5-10 ml) and a cannula.
3. Afterwards, take the products out of the cleaning bath and rinse at least three times thoroughly with water. With regard to products with lumens: rinse all instrument lumens five times at the beginning and or at the end of the contact time using a disposal syringe (minimum volume 5-10ml) and a cannula.
4. Check the products (see chapters on 'Checking' and 'Maintenance').

Disinfection:

5. Place the disassembled, cleaned and checked products in the disinfecting bath for the prescribed contact time and ensure that the products are sufficiently submerged. Take care that the products do not touch. With regard to products with lumens: rinse all instrument lumens five times at the beginning and or at the end of the contact time using a disposal syringe (minimum volume 5-10 ml) and a cannula.
6. Afterwards, take the products out of the disinfecting bath and rinse them thoroughly with water at least three times. With regard to products with lumens: rinse all instrument lumens five times at the beginning and or at the end of the contact time using a disposal syringe (minimum volume 5-10 ml) and a cannula.
7. Dry the products by blowing them dry using filtered pressurised air.
8. Wrap the products, if possible, straight after removal (see chapter on 'Packaging', if necessary, after they have been dried off completely in a clean place).

Proof of the general suitability for effective manual cleaning and disinfecting has been provided by an independent certified test laboratory using the cleaning agent Gigazyme and the disinfectant Gigasept Instru AF (Schulke & Mayr GmbH, Norderstedt). As part of the testing, the laboratory used the above described procedure.

Checking

After all products have been cleaned and/or cleaned/disinfected, check them for corrosion, damaged surfaces/bare patches, broken/chipped-off edges, deformations (e.g. bent rather than round) and impurities and eliminate damaged products (limited numbers for re-use see chapter on 'Re-use'). Products that are still contaminated need to be cleaned and disinfected once more.

Maintenance

Re-assemble disassembled products (see specific instructions). Instrument oils must not be used.

Packaging

Arrange the cleaned and disinfected products in the dedicated bur block/sterilisation tray. Wrap the products and/or the bur blocks/sterilisation trays using disposable sterilisation packaging (disposable or double packaging) and/or sterilisation containers that meet the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilisation (thermal stability up to 137 °C (279 °F) sufficient steam permeability)
- Sufficient protection of the products and/or sterilisation packaging against mechanical damage
- Regular maintenance according to the manufacturer's specifications (sterilisation container)
- Individual packaging: the packaging must be sufficiently large to ensure that the sealing is tension-free.

Sterilization

We only recommend the use of the sterilisation procedures listed below. Steam sterilisation

- Use of a fractional vacuum process or a gravitation process¹ (with sufficient product drying)
 - Steam sterilizer according to DIN EN 13060 and/or DIN EN 285
 - Validated according to DIN EN ISO 17665 (up to now: DIN EN 554/ ANSI AAMI ISO 11134) (valid IQ/OQ (commissioning) and product-specific performance assessment (FQ))
 - Maximum sterilisation temperature 134 °C (273 °F; plus tolerance according to DIN EN ISO 17665 (up to now: DIN EN 554/ANSI AAMI ISO 11134))
 - Sterilisation time (exposure time at the sterilisation temperature) at least 20 min at 121 °C (250 °F) and/or at least 52 min. at 132 °C (270 °F)/134 °C (273 °F)
- ¹Use of the less effective gravitation procedure is only admissible if the fractional vacuum process is unavailable.
² or 18 min (prion deactivation)

Proof of the general suitability of the products for effective steam sterilisation has been provided by an independent certified test laboratory using the steam sterilizer Systec V-150 (Systec GmbH Labor-Systemtechnik, Wetztenberg) as well as the fractional process and the gravitation process. During the testing, the laboratory took into account the typical hospital and practice conditions as well as the procedure described above.

As a rule, the flash sterilisation procedure is not admissible. Furthermore, do not use hot air sterilisation, radiation sterilisation, formaldehyde or ethylene oxide sterilisation or plasma sterilisation. In order to avoid stains and corrosion, the steam must be substance-free (see limit values included in DIN EN 13060). When sterilising several devices, the maximum load of the sterilising apparatus must not be exceeded (observe the manufacturer's instructions).

Storage

Prior to the first use of the device, the product should be stored in its original packaging at room temperature in dust- and humidity-free conditions. Subsequently, the products should be stored in appropriate hygienically maintained containers (protected from dust, humidity and recontamination). After sterilisation, the products need to be stored in sterilisation wrapping in a dry and dust-free place. Please note the shelf-life resulting from the validation of the sterilisation wrapping.

Material resistance

When choosing the cleaning and disinfecting agents ensure that they do not contain the following components:

- organic, mineral and oxidising acids (minimum admissible pH value 5.5)
 - strong alkaline solutions (maximum admissible pH value 8.5, neutral/enzymatic cleaning agent recommended)
 - organic solutions (e.g. alcohols, ether, ketones, benzines)
 - oxidants (e.g. hydrogen peroxides) - halogens (chlorine, iodine, bromine)
 - aromatic/halogenated hydrocarbons
- Never clean products, bur blocks and sterilisation trays using metal brushes or steel wool. Products, bur blocks and sterilisation trays should never be exposed to temperatures exceeding 137 °C (279 °F)!

Re-use

When it comes to products that are re-used, the end of their useful life tends to depend on their wear and tear as well as the use-related damage. Therefore, it is necessary to check the instrument after or prior to each re-use (see section on 'Checking'). Independently of this, it is the sole responsibility of the physician using the instrument to decide upon re-use based on the respective case as well as the potential wear and tear of the instrument. In this context, known restrictions as to the frequency of instrument use need to be observed.

If you fail to observe them, we do not accept any liability.

Special instructions

Drill extension:

- Disassemble instrument completely prior to cleaning and disinfecting
- Reassemble instrument prior to wrapping and sterilisation

Diamond products and ceramic grinding tools:

- Use particular care during the cleaning of the grinding surfaces and ensure that all residues are removed

Products equipped with a conduit to supply a cooling medium and other products with lumens (canals, drillings etc.):

- Actively rinse the lumen during pre-cleaning and manual cleaning and disinfection process

Supports:

- Disassemble instrument completely prior to cleaning and disinfecting
- Wrapping and sterilisation only in disassembled state.

Bur blocks/instrument trays:

- Cleaning and disinfecting only without products being loaded (products must not be cleaned and disinfected whilst they are in the bur block/instrument tray)
- Remove attachments prior to cleaning and disinfecting and clean and disinfect in disassembled state

The above instructions have been validated by the manufacturer of the medical devices who found them to be SUITABLE for preparing a medical device for re-use. It is up to the person in charge of the reprocessing to ensure that, based on the use of the correct equipment, material and personnel in the reprocessing facility, the actual reconditioning process produces the desired results. Normally, this requires the validation and routine monitoring of the procedure. Equally, each deviation from the instructions provided should be carefully checked for effectiveness and potential adverse consequences by the person in charge of reprocessing.

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**General application and safety instructions for the
medical device produced by MEISINGER
Issued: August 2013**

- MEISINGER medical products (podiatry) must only be used by podiatrists and/or the respective experts who, due to their training and experience, are intensely familiar with the use of these products and who have the corresponding expertise in the respective specialist fields.
- It is the sole responsibility of the operator in charge who, depending on the respective situation (indication), decides on the actual use of the MEISINGER products for each individual case.
- On request, training events are offered regarding the fields of application of MEISINGER foot care products. You will find further information on MEISINGER products on the internet at www.meisinger.de.
- All MEISINGER products have been developed for specific applications. Therefore, inappropriate use can result in the premature wear and tear of the products and put patients and users at risk.

Application

- In order to avoid damaging the instruments, they must be removed from the blister pack individually.
- It is essential to only use turbines as well as hand and angle pieces that are technically and hygienically faultless, maintained and cleaned.
- Rotating instruments need to be clamped as far down as possible with their speed set before applying them on the object. are used with the rotary instruments.
- Using the instruments for canting or leveraging should be avoided as it increases the risk of breakage.
- Depending on the application, it is recommended to use protective goggles as well as mouth protection while using the instruments.
- Users of diamond disks should use a disk protection device.
- Inappropriate use of the products leads to badly executed work and increased risk.
- A suction system should be used due to the development of dust when working with cutting and filing tools.
- In particular, users of hand tools should take care to use them gently and with consideration.
- The user must at all times avoid touching the instruments and parts without protection (protective gloves should be worn).

Use of pressure

- Users of the instruments should at all times avoid applying excessive pressure. This can damage the working part of the instruments and cause the cutting edges to break off. At the same time, it generates excessive heat.
- The use of excessive pressure when using grinding tools can cause the abrasive particles to break off or the instrument to become clogged and lead to heat generation.
- During polishing, excess pressure can lead to heat generation.
- Due to overheating, excess pressure can damage the texture or, due to broken off cutting edges, it can result in undesired effects. In such cases, even instrument breakage cannot be excluded.

Cooling

- In order to avoid excessive heat generation during preparation, a sterile-water/sodium chloride solution supplied via a permanent external feeding device should be used to ensure sufficient cooling during use of the instruments.

Storage, disinfection, cleaning and sterilization

- Unless there is explicit information to the contrary, all MEISINGER products are supplied in non-sterile packaging and, depending on the application, they need to be sterilised prior to use. Prior to their first use on the patient and immediately after each use, all products need to be disinfected and sterilised. Inappropriate cleaning and sterilising of the instruments can result in the patient being infected with harmful bacteria.
- You will find detailed instructions for the disinfecting, cleaning and sterilising of products in the Instructions for the reprocessing of medical devices produced by Hager & Meisinger GmbH. We would be happy to provide you with these instructions at your request. They are also available on the internet at www.meisinger.de.
- The products should be stored in appropriate, hygienically maintained containers. The same applies to sterilised instruments. The stored products must be protected from dust, humidity and recontamination. Instructions as to maximum storage duration must be adhered to.

Speed recommendations for rotary instruments

- Following the instrument-specific speed recommendations produces the best results.
- Exceeding the maximum admissible speed (rpm) when using long and pointed instruments tends to produce vibrations that can lead to the destruction of the instrument.
- When using working parts with diameters exceeding the thickness of the

shaft, excessive speed can release great centrifugal forces that may cause the shaft to bend and/or the instrument to break. Therefore, the maximum admissible rpm must never be exceeded.

- Please consult the manufacturer's information (see catalogue or www.meisinger.de) for the recommended and the maximum admissible speed ranges. Non-compliance with the maximum admissible speed puts safety at risk.
- Generally, the following rules apply:
– The larger the working part of an instrument the lower the speed
– The larger the working part of an instrument the greater the force of pressure

Discarding worn instruments and parts

- As a rule, MEISINGER medical products are re-usable unless there is explicit proof and information to the contrary. It is the sole responsibility of the operator using the devices to decide on the possibility of and the justification for multiple use as well as on the frequency of use, based on the respective case and the potential wear and tear of the products. If in doubt, it is always best to discard the products early and to replace them.
- Instruments are subject to wear and tear. This is particularly true for rotating tools. Therefore, our warranty for all rotating tools expires after their fifth use. If the tools are used more frequently, the manufacturer, Hager & Meisinger GmbH, can no longer guarantee a faultless cutting performance combined with minimum temperature change.
- Broken off cutting edges of instruments cause vibrations and great forces of pressure, which, in turn, leads to broken preparation corners and rough surfaces.
- Bare patches on diamond instruments indicate a lack of abrasive particles and can be a sign of blunt instruments. This leads to excessive temperatures during instrument use.
- Instruments that are bent and/or do not run true should be discarded forthwith.

Additional instructions regarding the use of trepanes

- When using trepanes, you have to proceed with particular care. For example, it is advisable not to exceed the recommended rpm speed ranges.
- In order to prepare for the actual use of a trepan, it should be set to produce counter-clockwise rotations creating a groove in the texture. Afterwards the trepan can be inserted into this groove and, using clockwise rotations, it can be moved further down.

Further comments:

- Due to statutory regulations, returned goods can, on principle, only be accepted if the complete batch number is provided. This number can be found on the product packaging.

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All rights reserved.
Reproduction, also by extract and reproducing photos, are only permitted with written authorization of Hager & Meisinger GmbH, Neuss.
All offers, orders, and deliveries are subject to Meisinger's "General Sales and Delivery Terms."
We reserve the right to modify our range of products and their design as well as to deviate from the illustrations and data shown.

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Certified since 1995 according to
DIN EN ISO 9001 and EN 46.000/13485